
## EUROPE ONLY

Distributed by:

AxoGen, Incorporated
13631 Progress Blvd., Suite 400
Alachua, Florida 32618
Phone: (386) 462-6880
Toll Free: (888) 296-4361
Fax: (386) 462-6801
E-Mail: CustomerCare@AxoGenInc.com
www.AxoGenInc.com

	<p>EC Representative:</p> <p>Cook Biotech Incorporated</p> <p>1425 Innovation Place</p> <p>West Lafayette, Indiana 47906 U.S.A.</p> <p>Phone: (765) 497-3355</p> <p>Toll Free: (888) 299-4224</p> <p>Fax: (765) 497-2361</p>
	<p>O'Halloran Road</p> <p>National Technical Park</p> <p>Limerick, IRELAND</p> <p>Phone: +353 61 334440</p>
	

©2013 AxoGen, Inc. AxoGuard® Nerve Protector and its logo are registered trademarks of AxoGen, Inc. U.S. Patents: 6,206,931; 7,652,077; 6,241,981; 6,358,284

Description
The AxoGuard® Nerve Protector is an implant that provides non-constricting protection for peripheral nerves. AxoGuard® Nerve Protector is designed to be an interface between the nerve and the surrounding tissue. AxoGuard® Nerve Protector is an extracellular matrix (ECM) and is fully remodelled during the healing process. When hydrated, AxoGuard® Nerve Protector is easy to handle, soft, pliable, non-frangible, and porous. AxoGuard® Nerve Protector is flexible to accommodate movement of the joint and associated tendons, and has sufficient mechanical strength to hold sutures. AxoGuard® Nerve Protector is provided sterile, for single use only, and in a variety of sizes to meet the surgeon's needs.

Indications for Use
The AxoGuard® Nerve Protector is intended for the repair of peripheral nerve discontinuities where gap closure can be achieved by flexion of the extremity. The AxoGuard® Nerve Protector is supplied sterile and is intended for single use.

**Rx ONLY** This symbol means the following:
**CAUTION: FEDERAL (U.S.A.) law restricts this device to sale by or on the order of a physician.**
This product is intended for use by trained medical professionals.

Contraindications
• The AxoGuard® Nerve Protector is derived from a porcine source and should not be used for patients with known sensitivity to porcine material.

**NOTE: This device not intended for use in vascular applications.**

**Precautions**
• This device is designed for single use only. Attempts to re-process, sterilize, and/or reuse may lead to device failure and/or transmission of disease.
• Do not resterilize device.
• Do not open or remove portions of device.
• Device is sterile provided the package is dry, unopened and undamaged. Do not use device if the package seal is damaged or open.

• Discard device if mishandling has caused possible damage or contamination, or if the device is past its expiration date.
• Do not suture device prior to rehydration.

**Potential Complications**

Possible complications can occur with any nerve repair surgical procedure including pain, infection, decreased or increased nerve sensitivity, and associated complications with anesthesia.

If any of the following conditions occur and cannot be resolved, careful removal of the device should be considered:
• Infection
• Allergic reaction
• Acute or chronic inflammation (initial application of surgical graft materials may be associated with transient, mild, localized inflammation)

**Storage**
The AxoGuard® Nerve Protector should be stored in a clean, dry location at room temperature.

**Sterilization**
This device has been sterilized with ethylene oxide.

**Suggested Instructions for Use**
**NOTE: These recommendations are designed to serve only as a general procedure. They are not intended to supersede institutional protocols or professional clinical judgment concerning patient care. Always handle AxoGuard® Nerve Protector using aseptic technique. Minimize contact with the device.**


1. Follow standard operating procedures for exposure and mobilization of the nerve. Determine the nerve diameter in millimeters (mm) using a suitable measuring instrument. If necessary, repair the nerve using standard opening procedures. Select an AxoGuard® Nerve Protector of sufficient diameter to account for normal edema following traumatic nerve injury. The AxoGuard® Nerve Protector diameter should be at least 1-2mm larger than the measured nerve diameter and long enough to cover the affected area.

2. Open the outer carton and remove the sterile pouch. Using standard aseptic technique, open the pouch and pass the inner tray to the sterile field for further handling.

3. Open the tray, and, if necessary trim the AxoGuard® Nerve Protector to the appropriate dimensions for covering the damaged portion of the nerve. The AxoGuard® Nerve Protector may be trimmed after hydration if desired.

4. Fill the pre-molded rehydration reservoir with room temperature sterile saline or sterile Lactated Ringer's solution. Hydrate the AxoGuard® Nerve Protector for 10 seconds or until the desired handling characteristics are achieved.

5. Position the AxoGuard® Nerve Protector around the nerve (Figure 1). If desired, gently flush the device with sterile saline or Lactated Ringer's solution to improve conformability with the nerve. Secure the device as necessary. This may include placing running sutures along the longitudinal slit to expose the nerve, and/or stay sutures through the nerve epineurium. See Figure 2 for completed repair.


6. Discard any unused portions of the AxoGuard® Nerve Protector according to institutional guidelines for biological waste. Do not resterilize.

**How Supplied**
AxoGuard® Nerve Protector is placed into a plastic tray and then inserted into a sterile pouch. The pouch is heat-sealed to provide a sterile barrier and has a peelable seal. Contents of the package are guaranteed sterile unless the package is opened or damaged. The AxoGuard® Nerve Protector and packaging do not contain natural rubber latex. Do not use AxoGuard® Nerve Protector if the peelable pouch appears to be open or damaged.

For additional information, to place an order or to report adverse events, contact:
AxoGen Customer Care: 888-AXOGEN (888-296-4361)
E-mail: CustomerCare@AxoGenInc.com

**Returned Goods Policy**
Authorization from AxoGen Customer Care must be obtained prior to returning device. Sterile product must be returned in unopened, undamaged cartons, packed to prevent damage.

**Symbols Used on Labeling**

	See instructions for use
	Expiration date
	Do not reuse after opening
	Lot number
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician
	Manufacturer
	Sterile unless package opened or damaged. Method of sterilization – ethylene oxide.

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