AxoGuard™
NerveConnector™

Distributed by:
AoxGen, Incorporated
13631 Progress Blvd., Suite 400
Cocke Bluffs 22015
Phone: (888) 462-6880
Toll Free: (888) 296-4361
CustomerCare@AoxGenInc.com
www.AxoGen.com

Manufactured by:
Cook Inland Incorporated
1425 Innovation Place
West Lafayette, Indiana 47906 U.S.A.
Phone: (765) 497-3355
Toll Free: (888) 299-4224
Fax: (765) 497-2361

EC Representative:
National Technological Park
Limerick, IR92
Phone: +353 61 334440
Fax: (353) 497-2361

©2013 AxoGen, Inc.
AxoGuard™ Nerve Connector and its logo are registered trademarks of AxoGen, Inc. U.S. Patents: 6,206,931; 7,652,077; 6,241,981; 6,358,284

Description

The AxoGuard™ Nerve Connector is a surgical implant designed to align and connect peripheral nerves. AxoGuard™ Nerve Connector is designed to be an interface between the nerve and the surrounding tissue. AxoGuard™ Nerve Connector is an extracellular matrix (ECM) and is fully remodeled during the healing process. When hydrated, AxoGuard™ Nerve Connector is easy to handle, soft, pliable, nonflexible, and porous. AxoGuard™ Nerve Connector is flexible to accommodate movement of the joint and associated tendons, and has sufficient mechanical strength to hold sutures. AxoGuard™ Nerve Connector is provided sterile, for single use only, and in a variety of sizes to meet the surgeon's needs.

Indications for Use

The AxoGuard™ Nerve Connector is intended for the repair of peripheral nerve discontinuities where gap closure can be achieved by flexible alignment. AxoGuard™ Nerve Connector is supplied sterile and is intended for single use.

[Rx ONLY] This symbol means that you should be cautious. CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician. This product is intended for use by trained medical professionals.

Contraindications

The AxoGuard™ Nerve Connector is derived from a porcine source and should not be used in patients with known sensitivity to porcine material.

NOTE: This device is not intended for use in vascular applications.

Precautions

- This device is designed for single use only. Attempts to reprocess, sterilize, and/or reuse may lead to device failure and/or transmission of disease.
- Do not sterilize device.
- Discard all open and unused portions of device.
- Device is sterile provided the package is dry, unopened and undamaged. Do not use if the package seal is damaged or if the device is damaged.
- Discard device if mishandling has caused possible damage or contamination, or if the device is past its expiration date.
- Do not reuse device prior to rehydration.

Potential Complications

Possible complications can occur with any nerve repair surgical procedure including pain, infection, decreased or increased nerve sensitivity, and complications associated with use of anesthesia.

• Allergic reaction

• Acute or chronic inflammation (initial application of surgical graft materials may be associated with transient, mild, localized inflammation)

Storage

AxoGuard™ Nerve Connector should be stored in a clean, dry location at room temperature.

This device has been sterilized with ethylene oxide.

Suggested Instructions for Use

NOTE: These recommendations are designed to serve only as a general procedure. They are not intended to supersede institutional protocols or professional clinical judgment concerning patient care. Always handle AxoGuard™ Nerve Connector using aseptic technique. Minimize contact with latex gloves.

1. Follow standard operating procedures for exposure and mobilization of the nerve (see Figure 1). Determine the nerve diameter in millimeters (mm) using a suitable measuring instrument. Select an AxoGuard™ Nerve Connector of sufficient diameter to account for normal edema following traumatic nerve injury and to allow easy insertion of the nerve stumps into the device.

 **Figure 1**

2. Open the outer carton and remove the sterile pouch. Using standard aseptic technique, open the pouch and pass the inner tray to the sterile field for further handling.

3. Open the tray and, if necessary, trim the AxoGuard™ Nerve Connector to a length that is appropriate for the nerve gap. Allow at least 2mm on the proximal and distal stumps for insertion into the connector. The AxoGuard™ Nerve Connector may be trimmed after hydration if desired.

4. Fill the pre-moiled rehydration reservoir with room temperature sterile saline or sterile Lactated Ringer's solution. Hydrate the AxoGuard™ Nerve Connector for 10 seconds or until the desired hydration characteristics are achieved.

5. Hemostasis of both nerve stumps must be achieved prior to beginning the entubulation procedure. In the event a tourniquet is used, release the tourniquet and achieve hemostasis before beginning the entubulation procedure.

6. Use non-absorbable suture to secure the AxoGuard™ Nerve Connector in place. Pass the suture through the wall of the tube from outside to inside at least 1mm from the edge (see Figure 2). Then pass the suture through the epineurium of one nerve stump at a distance of at least 1 mm from the cut nerve face. Reverse the suture and pass it through the wall of the AxoGuard™ Nerve Connector by pulling the suture so that the nerve stump is drawn into the tube. Tie the suture so that it is secure but avoid generating tension in the suture itself. If desired, add an additional suture to secure the device.

 **Figure 2**

7. Use a syringe to gently flush the lumen of the AxoGuard™ Nerve Connector with sterile saline or Lactated Ringer's solution (Figure 3). Repeat the flushing procedure described above with the opposite nerve stump (Figure 4). Fill the interior of the AxoGuard™ Nerve Connector with saline or Lactated Ringer's solution (Figure 5). See Figure 6 for a completed repair.

 **Figure 3**

 **Figure 4**

 **Figure 5**

 **Figure 6**

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