**DESCRIPTION**  
Avance Nerve Graft is a processed nerve allograft intended for bridging nerve discontinuities to support regeneration. It is decontaminated and cleaned extracellular matrix from donated human peripheral nerve. The cleaning process preserves the inherent microtubular structure and associated cytoskeletal proteins, including microtubules, cellular debris and certain proteins such as chondroitin sulfate proteoglycans (CSPG).

Avance Nerve Graft is implanted and connects the proximal and distal ends to a transected nerve. The processed nerve allograft serves as a scaffold allowing regenerating axons to grow into the patient's distal nerve tissue toward the target muscle or skin.

Avance Nerve Graft offers structural characteristics and handling similar to nerve allograft; pliability of soft tissue, an intact epineurium to suture the nerve graft in place, and intact endoneurial tubes for the axons to grow through.

**POTENTIAL COMPLICATIONS**  
The following adverse events are anticipated and may potentially occur during treatment:  
• mild incisional redness,  
• tenderness of surgical area,  
• mild edema of surgical area,  
• controllable pain at surgical area,  
• decreased sensation,  
• numbness.

**REGULATORY CLASSIFICATION**  
Avance Nerve Graft is processed and distributed in accordance with US FDA requirements for Human Cellular and Tissue-based Products (HCT/P) under 21 CFR Part 1271 regulations, US State regulations of the American Association of Tissue Banks (AATB). Additionally, international regulations are followed as appropriate.

Avance Nerve Graft is to be dispensed only by or on the order of a licensed physician.

**INDICATIONS FOR USE**  
Avance Nerve Graft is processed nerve allograft (human) intended for the surgical repair of peripheral nerve discontinuities to support regeneration across the defect.

**CONTRAINDICATIONS**  
Avance Nerve Graft is contraindicated for use in any patient in whom soft tissue implants are contraindicated. This includes any pathology that would limit the blood supply and compromise healing or evidence of a current infection.

**WARNINGS**  
Careful donor screening, laboratory testing, tissue processing, and gamma irradiation have been utilized to minimize the risks of transmission of relevant communicable diseases to the patient. As with any processed human donor tissue, Avance Nerve Graft cannot be guaranteed to be free of all pathogens.

Do not use or resterilize Avance Nerve Graft and do not refreeze the graft once it has been thawed.

**DONOR RECOVERY AND SCREENING**  
After consent for donation is obtained, surgical recovery of the peripheral nerve tissue is performed in FDA registered and state licensed facilities. Avance Nerve Graft is processed in FDA registered and state licensed facilities. Avance Nerve Graft is processed in FDA registered and state licensed facilities. Avance Nerve Graft is processed in FDA registered and state licensed facilities. Avance Nerve Graft is processed in FDA registered and state licensed facilities.

**Single European Code (SEC) for Tissues**  
If applicable SEC Code is found on or attached to the outer carton.

Each donor is tested and shown to be negative or nonreactive for the following:  
• Human Immunodeficiency Virus (HIV) Type 1 Antibody  
• Human Immunodeficiency Virus (HIV) Type 2 Antibody  
• Hepatitis C Virus (HCV) Antibody  
• Hepatitis B Virus (HBV) Surface Antigen  
• Hepatitis B Virus (HBV) Core Antibody (total)  
• Syphilis Rapid Plasma Reagin or Treponemal Specific Assay  
• Human Immunodeficiency Virus (HIV) Nucleic Acid Test (NAT)  
• Hepatitis C Virus (HCV) Nucleic Acid Test (NAT)  
• Hepatitis B Virus (HBV) Nucleic Acid Test (NAT)

**REFERENCES**  
1. ISO 11317-1:2006/Am2 2:2018 Sterilization of health care products - Radiation guidelines  
2. ISO 10993:2003 Biological evaluation of medical devices

**INQUIRIES**  
For additional information, to place an order, or to report errors, accidents or adverse reactions, contact:  
Axogen Customer Care  
Phone: 888.496.9111 (888.296.4361) or 386.462.8800  
Email: [customer-care@axogeninc.com](mailto:customer-care@axogeninc.com)

Customers outside of US:  
Contact Axogen authorized distributor servicing your facility, email [Axogen Customer Care](mailto:Axogen Customer Care) or contact Axogen directly in the US at 386.462.8800

Symbols used on Product Packaging:

REF Catalogue number  
LOT Lot number  
Temperature limit  
STERILE R Sterilized using irradiation  
Do not resterilize  
Do not reuse  
Consult instructions for use  
Use-by date  
D Diameter  
L Length

The Medical Director of Axogen (US state licensed) has determined that the tissue is suitable for transplantation in humans. Records of all testing and medical releases are maintained by Axogen.

**PROCESSING**  
Avance Nerve Graft is processed in controlled environments using Good Tissue Practice (GTP) methods designed to prevent contamination and cross-contamination of tissue. Processing involves the use of proprietary physiological buffers, enzyme and surfactants and the processed tissue may contain traces of these processing agents. The cleaning process preserves the inherent microtubular structure of the native nerve extracellular matrix while cleaning the cells, cellular debris and certain proteins such as chondroitin sulfate proteoglycans (CSPG). After cleaning, the tissue is decontaminated and sterilized using gamma irradiation in accordance with ISO 11317 guidelines.

Avance Nerve Graft has been tested in accordance with ISO 10993 standards. The test results demonstrated that the processed nerve allograft is biocompatible.

The processed nerve allograft tissue will naturally vary in color from white, off-white, pink, pale pink, and yellow to pale yellow. Occasional dark spots or localized discoloration is a normal occurrence.

**HOW SUPPLIED**  
Avance Nerve Graft is placed into a plastic tray and then inserted into chevron pouches. Each chevron pouch is heat-sealed to provide a sterile barrier and each pouch has a chevron seal. The outer chevron pouch is foil to provide a moisture barrier. Approximate graft lengths and diameters are listed on the package label. Avance Nerve Graft is irradiated and supplied frozen. Contents of the foil package are sterile unless the package is opened or damaged.

**TRANSPORT AND STORAGE**  
Avance Nerve Graft is shipped frozen on dry ice via validated shipping containers. Upon receipt, verify the contents, verify the temperature of the shipment and place the grafts in a cooler or below -40°C (-40°F). Keep product frozen until use. If no dry ice is present, contact Axogen Customer Care as product may have thawed. Expiration date is three (3) years from date of packaging, provided that the Avance Nerve Graft has been stored at temperatures at or below -40°C (-40°F).

**See product label for expiration date. Expiration date of the product is the end of the month for the month and year outlined on the product label.**

**Temporary storage conditions**  
Avance Nerve Graft may be stored at -20°C to -40°C (-4°F to -40°F) for up to 3 months. If stored at -20°C (-4°F) for longer than 3 months, the grafts should be stored at this temperature range must then be transferred to -40°C (-40°F) or colder freezer, used or discarded.

**It is the responsibility of the health care institution to track the expiration date of the Avance Nerve Graft and ensure that the product is stored properly.**

**WARNING**  
If the outer foil chevron pouch and/or inner Tyevek® chevron pouch is compromised (shows evidence of being torn or opened in any manner), DO NOT USE the Avance Nerve Graft and notify Axogen Customer Care immediately.

**INSTRUCTIONS FOR USE**  
1. Follow standard operating procedures for exposure and mobilization of the injured peripheral nerve.  
2. Determine the injured nerve diameter and length in millimeters (mm) using a sutureless measuring instrument.  
3. Select Avance Nerve Graft(s) of comparable diameter to match the native nerve and of sufficient length to ensure a tension free repair.  
4. To prepare the Avance Nerve Graft:  
a. Remove the foil chevron pouch containing Avance Nerve Graft. Instructions for Use, patient record labels and Tissue Utilization Report (TUR) from the package.  
b. Compare the distinct lot number on the foil chevron pouch with the lot number on the package. If the numbers do not match, DO NOT USE the product and notify Axogen Customer Care immediately.  
c. Using standard aseptic technique, peel open the outer foil chevron pouch and pass the inner Tyevek® chevron pouch to the sterile field for further handling.  
d. Open the Tyevek® chevron pouch and remove the plastic tray.  
e. Open the plastic tray and fill the pre-sterilized thawing reservoir with room temperature sterile saline or sterile Lactated Ringers' Solution (LRS). Do not heat the graft or add heated saline or LRS to the graft.  
f. Allow Avance Nerve Graft to thaw completely before use which will take approximately 5 – 10 minutes. Once thawed, Avance Nerve Graft is soft and pliable throughout. Avance Nerve Graft must be kept at room temperature for a maximum of 12 hours. NEVER IMPLANT A PARTIALLY OR FULLY FROZEN PRODUCT.  
g. Handle Avance Nerve Graft by outer most epineurium and avoid crimping or crushing the nerve.  
h. Implant Avance Nerve Graft using the same tensionless surgical technique used when implanting a nerve allograft. Either end of the processed nerve allograft can be coapted to the proximal stump of the host nerve.  
i. Destroy any thawed allograft tissue not used in the surgical procedure in accordance with local, state and federal regulatory requirements for disposal of human tissue.  
j. Complete and send the Tissue Utilization Report (TUR) back to Axogen.

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