Distributed by:  
 Axogen Corporation  
 13631 Progress Blvd., Suite 400  
 Alachua, Florida 32015  
 Phone: (386) 462-6880  
 Toll Free: (888) 296-4361  
 Fax: (386) 462-6801  
 CustomerCare@axogeninc.com  
 www.Axogeninc.com

Manufactured by:  
 Cook Biotech Incorporated  
 1425 Innovation Place  
 West Lafayette, Indiana 47906 U.S.A.

AxoGuard™ Nerve Connector and its logo are trademarks of Axogen Corporation, U.S. Patent Pending.

**Description**  
 The AxoGuard Nerve Connector is a surgical implant designed to align and connect peripheral nerves. AxoGuard Nerve Connector is designed to be an interface between the nerve and the surrounding tissue, thus providing a protective environment for the growing nerve healing. AxoGuard Nerve Connector is comprised of an extracellular matrix (ECM) and is fully remodeled during the healing process. Animal model data suggests that the AxoGuard Nerve Connector is fully vascularized in situ and gradually remodeled into a connective tissue structure similar to nerve epineurium over several months. When hydrated, AxoGuard Nerve Connector is easy to handle and manipulate, and porous. AxoGuard Nerve Connector is flexible to accommodate movement of the joint and associated tendons, and has sufficient mechanical strength to hold only, and AxoGuard Nerve Connector is provided sterile, for single use only, and in a variety of sizes to meet the surgeon's needs.

**Indications for Use**  
 The AxoGuard Nerve Connector is indicated for the repair of peripheral nerve discontinuities with gaps up to 5 mm. The AxoGuard Nerve Connector is supplied sterile and is intended for single use.

**[Rx ONLY]** This symbol means the following:  
**CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.**

This product is intended for use by trained medical professionals.

**Contraindications**  
 The AxoGuard Nerve Connector is derived from a porcine source and should not be used for patients with known sensitivity to porcine material.

**NOTE: This device is not intended for use in vascular applications.**

- Precautions**
- This device is designed for single use only. Attempts to reprocess, sterilize, and/or reuse may lead to device failure and/or transmission of disease.
  - Do not resterilize device.
  - Discard all open and unused portions of the device.
  - Devices is sterile provided the package is dry, unopened and undamaged. Do not use device if the package seal is damaged or open.
  - Discard device if mishandling has caused possible damage or contamination, or if the device is past its expiration date.
  - Do not store device prior to rehydration.

**Potential Complications**  
 Possible complications can occur with any nerve repair surgical procedure. These include, but are not limited to, wound healing, nerve sensitivity, and complications associated with use of anesthesia.

Any of the following conditions occur and cannot be resolved, careful removal of the device should be considered.

- Acute or chronic inflammation (initial application of surgical graft materials may be associated with transient, mild, localized inflammation)
- Allergic reaction
- Infection

**Storage**  
 The AxoGuard Nerve Connector should be stored in a clean, dry location at room temperature.

**Sterilization**  
 This device has been sterilized with ethylene oxide.

**Suggested Instructions for Use**  
**NOTE: Always handle AxoGuard Nerve Connector using aseptic technique. Minimize contact with latex gloves.**

1. Follow standard operating procedures for exposure and mobilization of the nerve (see Figure 1). Determine the nerve length in millimeters (mm) using a suitable measuring instrument. Select an AxoGuard Nerve Connector of sufficient diameter and length to allow for mobilization of the nerve and to allow easy insertion of the nerve stumps into the connector. The AxoGuard Nerve Connector may be trimmed after hydration if desired.

2. Open the outer carton and remove the sterile pouch. Using standard aseptic technique, open the pouch and pass the inner tray to the sterile field for further handling.

3. Open the tray and, if necessary, trim the AxoGuard Nerve Connector to a length that is appropriate for the nerve gap. Allow for an additional 5 mm nerve gap and at least 2 mm on the proximal and distal stumps for insertion into the connector. The AxoGuard Nerve Connector may be trimmed after hydration if desired.

4. Fill the pre-molded rehydration reservoir with room temperature sterile saline or sterile Lactated Ringer's solution. Hydrate the AxoGuard Nerve Connector for 10 seconds or until the desired handling characteristics are achieved, but not more than 20 minutes.

5. Hemostasis of both nerve stumps must be achieved prior to beginning the entubulation procedure. In the event a tourniquet was used, release the tourniquet and achieve hemostasis before beginning the entubulation procedure.

6. Use non-absorbable suture to secure the AxoGuard Nerve Connector in place. Pass the suture through the wall of the tube from outside to inside, at least 1 mm from the edge. See Figure 2. Then, pass the suture transversely through the epineurium of one nerve stump at a distance of at least 1 mm from the cut nerve site. Reverse the suture and pass it through the wall of the AxoGuard Nerve Connector by pulling the suture through the tube. Make sure the suture is drawn into the tube. Tie the suture so that it is secure but avoid generating tension in the suture itself. If desired, add an additional suture to secure the device.

7. Use a syringe to gently flush the lumen of the AxoGuard Nerve Connector with sterile saline or Lactated Ringer's solution (Figure 3). Repeat the suturing procedures described above with the opposite nerve stump (Figure 4). Fill the interior of the AxoGuard Nerve Connector with saline or Lactated Ringer's solution (Figure 5). See Figure 6 for complete repair.

**NOTE: Clinical literature has shown that better outcomes are found with using conduits with a gap of 5 mm or less.<sup>1-3</sup>**

**Figure 1**

**Figure 2**

**Figure 3**

**Figure 4**

**Figure 5**

**Figure 6**

**How Supplied**  
 AxoGuard Nerve Connector is provided in a plastic tray and outer sterile pouch. The pouch is heat-sealed to provide a sterile field and has a peelable seal. Contents of the package are guaranteed sterile unless the package is opened or damaged. The AxoGuard Nerve Connector and packaging do not contain natural rubber latex. Do not use if the AxoGuard Nerve Connector pouch appears to be open or damaged.

**References**  
 1. Kokkalis ZT, Puc S, Small GA, Weiser RW, Venouzou AI, Soteranos DG. Assessment of processed porcine extracellular matrix as a protective barrier in a rabbit nerve repair model. *J Reconstr Microsurg*. 2011;27(1):19-28.

2. Isaacs J, Safa B, Evans PJ, Greenberg J. Technical assessment of connector assisted nerve repair. *J Hand Surg*. 2016;41(7):760-766.

3. Duccic I, Safa B, DeVlinny E. Refinements of nerve repair with connector-assisted coaptation. *Microsurgery*. 2017;37(3):256-263. Inquiries

For additional information, to place an order, or to report adverse event contact:  
**Axogen Customer Care: 888-AXOGEN (888-296-4361)**

**Returns and Policy**  
 Authorization from Axogen Customer Care must be obtained prior to returning product. Sterile product must be returned in unopened, undamaged cartons, packed to prevent damage.

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**Do not reuse**

**Batch Code**

**Caution: Federal (USA) law restricts this device to sale by or on the order of a physician**

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